


GOLDEN RF MICRONEEDLING DEVICE INFORMATION



NEEDLES



Contacts

Needle depth range

Treatment area



9 pins

0.5-4.0mm

Around the eyes



20 pins

0.5-6.0mm

Armpit



25 pins

0.5-4.0mm

Neck




49 pins

0.5-3.0mm

Face
Abdomen

Each probe has 1200 rounds, and the number is 0 when used up.
The probes are replaced one after another and cannot be used interchangeably.

HOW IT WORKS



- 16, 25 and 49 needle array
- Adjustable depth 0.5mm - 4.0mm
- Used for treating more tissue in a single treatment, without sacrificing the epidermis

TREATMENTS BEFORE AND AFTER

ACNE



BEFORE

AFTER

WRINKLES



BEFORE

AFTER

NECK WRINKLES



BEFORE

AFTER

STRETCH MARKS





















BEFORE

AFTER

- Improve oil balance
- Improves large pores
- Improve acne
- Improve different degrees of acne pits

- Improves facial sagging and wrinkles
- Improve traumatic scars
- Improve stretch marks
- Improve fundus elasticity loss and relaxation

DEVICE COMPARISON

	RF GOLD	DERMAROLLER	DERMAPEN
			
Needle Output Automatically			
Adjustable Depth			
Optional Treatment Head			
Higher Efficiency			
RF Energy Output			
More Effective	